

SALIVA ALCOHOL TEST STRIP

Mission® Saliva Alcohol Test Strips, a rapid saliva test for the determination of relative Blood Alcohol Concentration (BAC). Results in 2 mins.

BROAD AND COMPLETE DETECTION

- Single strip detects multiple concentration of 0.02-0.30% BAC
- Unique strip with all-in-one detection offers cost effective testing of a single sample
- 18 months shelf life

FASY TO USE

- Non-invasive specimen collection allows for on-site testing
- No additional equipment or training necessary
- Easy to interpret results in 2 minutes.

ACCURATE AND RELIABLE

- Quick and easy detection
- Suitable for use in homes, school, workplaces, rehabilitation centers, hospitals, military, law enforcement units and other government agencies.

CONVENIENT

- Cost-effective method compared to other alcohol detection test
- Compact and individually wrapped packaging for use anytime and anywhere.



PARAMETERS	SPECIFICATION
B l ow Bag	N/A
Technology	Enzymatic Reaction
Specimen	Saliva
Blood Alcohol Concentration (BAC)	0.02%-0.30%
Measuring Units	%BAC
Reading Time	2 Minutes
Reading Stability	3 Minutes
Storage Temperature	2-27°C (36-81°F)
Shelf Life	18 Months
Package Size (Test)	25



