

Baso

HEMATOLOGY STAIN



DIFF QUICK STAIN

Diff-Quik Stain is mainly used for blood cell smear, bone marrow smear, vaginal discharge (gynecological vaginal discharge) smear, exfoliated cells smear staining. Available in 3 x 500ml/kit.



WRIGHT GIEMSA STAIN

This kit is for classification of hematocytes in blood and marrow. Available in 4 x 250 ml/kit.



RETICULOCTE STAIN

Reticulocyte Stain is intended for the intravital staining of reticulocytes and it is different from the traditional method which requires counterstain to distinguish the reticulocyte from mature erythrocytes. The background of this stain is bright and very clear, and is not affected by excessive staining. For an extended staining, the background will become brighter. Available in 6 x 20ml / kit.



IRON STAIN

This kit is for staining bone marrow cell and bloodcell smear.

The iron staining procedure utilizes the Prussian Blue stain for ferric iron to assess bone marrow iron stores. This procedure is particularly helpful when evaluating patients with anemia, iron overload and myelodysplasia. Available in 4 x 100ml / kit.



TruLaboratories
Corporation

+632 82-TRULAB (878522)
+63917.878.2288 • +63925.878.2288
www.trulaboratories.com